




## Titration head / 2 x SGJ 14

6.1458.070

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.

納品範囲 6.1458.070

| Qt.   | Order no.  | 説明                    |
|-------|------------|-----------------------|
| 4 PCS | 6.1808.000 | Connection bushing M6 |



3 PCS

6.2726.090

拡散防止バルブ

拡散防止用ビュレットチップ6.1541.XXX、 6.1543.130、  
6.1543.140、 6.1543.150、 6.1543.160、 6.1543.200



---

2 PCS

6.2740.020

Spray nozzle



# オプションアクセサリー

| Order no.  | 説明                                                                            |                                                                                     |
|------------|-------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.1829.050 | Replacement tubings for retracted tubings in titration heads and robotic arms |  |