



Robotic arm with holder for titration head and sensor for 786 Swing Head, right swinging

6.1462.170

オプションアクセサリ

Order no.	説明
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.



6.1458.020 Titration head / 4 x M10 connectors



6.1458.030 KF titration head / 4 x M10 connectors



6.1458.040 Titration Head, 3 x SGJ 14



6.1458.070 Titration head / 2 x SGJ 14

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.

