



サンプルラック 28 x 250 mL、(ポリ塩化ビニル製)

6.2041.820

28個のサンプルビーカーに適した 815 Robotic Sample Processor、855 Robotic Titrosampler「Basic」のためのサンプルラック。以下のサンプルビーカーの使用が可能です：
6.14332.320、6.1432.323、6.1453.220、6.1453.250。

合成樹脂: ポリ塩化ビニル (PVC)

オプションアクセサリ

Order no.	説明
6.1432.320	透明ガラス製サンプルビーカー 250 mL、1 個
	透明ガラス製サンプルビーカー 250 mL、サンプルラック 6.2041.310、6.2041.400、6.2041.430、6.2041.510、 6.2041.820、6.2068.020 および 6.2068.030 に適合、1 個



6.1453.220 プラスチック製 (PP) サンプルビーカー 200 mL、1 個

プラスチック製 (ポリプロピレン、PP) サンプルビーカー 200 mL、サンプルラック 6.2041.310、6.2041.400、6.2041.430、6.2041.510、6.2041.820、6.2068.020 および 6.2068.030 に適合、1 個



6.1453.250 プラスチック製 (PP) サンプルビーカー 250 mL、1 個

プラスチック製 (ポリプロピレン、PP) サンプルビーカー 250 mL、サンプルラック 6.2041.310、6.2041.400、6.2041.430、6.2041.510、6.2041.820、6.2068.020 および 6.2068.030 に適合、1 個



6.1458.010 Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.



6.1458.040 Titration Head, 3 x SGJ 14



6.1462.060 Robotic arm with holder for titration head, left swinging



6.1462.070 Robotic arm with holder for titration head, right swinging



6.1462.160 Robotic arm with holder for titration head and sensor for 786
Swing Head, left swinging



6.1462.170 Robotic arm with holder for titration head and sensor for 786
Swing Head, right swinging

