



Buffer solution pH 9.00 (30 x 30 mL)

6.2307.220

オプションアクセサリ

Order no.	説明
6.2307.230	Buffer solutions pH 4, 7 and 9



6.2307.200 Buffer solution pH 4.00 (30 x 30 mL)



6.2307.210 Buffer solution pH 7.00 (30 x 30 mL)



6.2307.100 Buffer solution pH 4 (500 mL)



6.2307.110 Buffer solution pH 7 (500 mL)



6.2307.120 Buffer solution pH 9 (500 mL)



6.2325.000 pHit kit

