



pHit kit

6.2325.000

納品範囲 6.2325.000

Qt.	Order no.	説明
2 PCS	6.2008.040	保存容器

ホルダー 6.2008.050 とともに使用すると 807 Dosing Unit
の電極用のホルダーになります。



オプションアクセサリ

Order no.	説明
6.2615.050	Electrode stand



6.2307.230 Buffer solutions pH 4, 7 and 9



6.2307.200 Buffer solution pH 4.00 (30 x 30 mL)



6.2307.210 Buffer solution pH 7.00 (30 x 30 mL)



6.2307.220 Buffer solution pH 9.00 (30 x 30 mL)



6.2307.100 Buffer solution pH 4 (500 mL)



6.2307.110 Buffer solution pH 7 (500 mL)



6.2307.120 Buffer solution pH 9 (500 mL)

