



Equipment for pH Adjustment at the 845 Eluent Synthesizer

6.5329.000

納品範囲 6.5329.000

Qt.	Order no.	説明
1 PCS	6.1446.000	SGJ stopper / B-14/(15)



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1602.190


Eluent Bottle Head, 1 x SGJ 14, 4 x M6



1 PCS

6.1838.010

Connecting Tubing to 845, 85 cm



1 PCS

6.2148.010

Remote Box MSB



1 PCS

6.2726.090

拡散防止バルブ

拡散防止用ビュレットチップ6.1541.XXX、6.1543.130、
6.1543.140、6.1543.150、6.1543.160、6.1543.200



1 PCS

6.3026.210

Exchange Unit 10 mL

