




# IC equipment for Liquid Handling with the Dosino

6.5330.040

納品範囲 6.5330.040

| Qt.   | Order no.  | 説明                       |
|-------|------------|--------------------------|
| 1 PCS | 6.1602.160 | 溶離液および試薬のためのボトルキャップGL 45 |

溶離液ボトル、試薬ボトルおよび洗浄ボトル用、乾燥管および吸引チューブ用コネクタ付き。



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



---

1 PCS

6.1618.020

Thread adapter S 40 to GL 45



---

1 PCS

6.1619.000

807 Dosing Unit 用乾燥管



---

1 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube



1 PCS

6.1805.120

FEP tubing / M6 / 100 cm



1 PCS

6.1819.100

FEP aspiration tubing to 2 L bottle



1 PCS

6.1831.030


PEEK capillary 0.75 mm i.d. / 3 m



1 PCS

6.2057.210

Holder for Dosino to IC instruments



1 PCS

6.2744.010

Pressure screw 5x



---

2 PCS

6.2744.080

M6 thread / UNF 10/32 adapter

