



ICアクセサリー: 溶離剤生成モジュールの追加 溶離剤

6.5330.090

溶離液生成モジュールの付属品セット 溶離剤拡張のための全てのパーツを含みます。レベルセンサーは用途に応じてご注文ください。

納品範囲 6.5330.090

Qt.	Order no.	説明
1 PCS	6.1602.160	溶離液および試薬のためのボトルキャップGL 45
溶離液ボトル、試薬ボトルおよび洗浄ボトル用、乾燥管および吸引チューブ用コネクタ付き。		



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS

6.1618.020

Thread adapter S 40 to GL 45



1 PCS

6.1618.050

Thread adapter 40 mm to GL 45



1 PCS

6.1621.000

Container 10 L



1 PCS

6.1621.100

Cap for 10-L container



1 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube



1 PCS

6.1626.000

Bottle attachment GL 45 for level sensor



4 PCS

6.1805.530

FEP tubing / M6 / 2 m



1 PCS

6.1808.280

807 Dosing Unit ポート 4 アダプター、M6 内側

FEP チューブ / M6 との接続のため、807 Dosing Unit のポート 4 に使用するアダプター。



1 PCS

6.1815.010

Spiral band 0.5 m



1 PCS

6.1819.100

FEP aspiration tubing to 2 L bottle



1 PCS

6.1819.110


FEP aspiration tubing to canister



1 PCS

6.2057.210

Holder for Dosino to IC instruments



1 PCS

6.2151.060

Cable for Level Control and Eluent Production Module



1 PCS

6.2744.200

Adapter UNF 10/32 outer / M6 inner



1 PCS

6.1580.250

807 Dosing Unit 50 mL、付属品なし





オプションアクセサリ

Order no.	説明	
6.2769.000	FULL sensor for level sensor	
6.2769.100	EMPTY sensor for level sensor for 1 L bottles	
6.2769.110	EMPTY sensor for level sensor for 2 L bottles	
6.2769.120	EMPTY sensor for level sensor for 10 L canisters	