



IC機器: インライン透析 Low Volume

6.5330.200

高速インライン透析のための付属品セット。858 Professional Sample Processorおよび追加の2チャンネルタイプのペリスタリックポンプとの使用に。

納品範囲 6.5330.200

Qt.	Order no.	説明
2 PCS	6.2744.070	Pressure screw short



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS

6.1831.160

PEEK capillary, 0.5 mm ID / 70 cm



2 PCS

6.2744.034

カップリング チューブニップル/UNF 10/32、2個


圧力調整ネジとポンプチューブを接続。2個。ペリスタリックポンプ付きIC装置用。



4 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device



1 PCS

6.1602.150

Bottle attachment / GL 45 - 3 x UNF 10/32



8 PCS

6.1803.040

PTFE capillary 0.5 mm i.d. / 1 m



2 PCS

6.1826.340

Pump tubing LFL (black/black), 3 stoppers



1 PCS

6.2755.000

Tubing cartridge



1 PCS

6.1831.120

PEEK capillary 0.25 mm i.d., 45 cm



1 PCS

6.2744.000

Pressure screw PVDF



2 PCS

6.1826.320

Pump tubing LFL (orange/yellow) , 3-stopper



1 PCS

6.2831.010


Plastic tweezers



1 PCS

6.2057.130

Dialysis cell holder



1 PCS

6.2714.010

Dialysis membrane (cellulose acetate)



1 PCS

6.2729.200

透析セル Low Volume

4 個のストッパー 6.2744.060 を含む。高速インライン透析用。

