



IQ/OQ Kit for IC

6.5333.000

納品範囲 6.5333.000

Qt.	Order no.	説明
1 PCS	6.1803.050	PTFE capillary 0.5 mm i.d. / 20 cm



1 PCS 6.1831.170 OQ capillary



1 PCS 6.2322.030 IC のための稼働性能適格性確認 (OQ) 標準液

IC (6.5333.000) 用IQ/OQキットの一部。個別に注文することはできません。



1 PCS 6.2744.040 2 x UNF 10/32 coupling



1 PCS 6.2744.060 Threaded stopper



