



Transfer equipment, 0.5 - 2 mL, for 789/815

6.5619.000

納品範囲 6.5619.000

Qt.	Order no.	説明
1 PCS	6.1446.180	Stopper, PTFE, SGJ 29, for transfer set



1 PCS

6.1448.000

Septum stopper B-14/15



2 PCS

6.1448.020

Septum 16 mm, 5 pieces



1 PCS

6.1462.090

Robotic arm with Luer-Lock adapter for 786
Swing Head, right swinging



1 PCS

6.1464.320

KF titration vessel / 80-250 mL / coulometric



1 PCS

6.1562.110

Pipetting tubing / 2 mL



1 PCS

6.1608.090


Bottle / 250 mL / GL 45



1 PCS

6.2001.100

Base plate, transfer set



1 PCS

6.2013.010

Clamping ring



1 PCS

6.2701.040

GL 18 screw cap with hole



1 PCS

6.2751.120

Safety shield Transfer set to 789/815



1 PCS

6.2816.100

Needle with vent



1 PCS

6.2055.110

ボトル2本 / Dosino のためのボトルホルダー、もしくはDosinoホルダー

250 mL のボトル2本もしくは Dosino 2つのためのホルダー

