




Consumable Kit Titration

6.5702.010

納品範囲 6.5702.010

Qt.	Order no.	説明
1 PCS	6.1543.200	Titration tip / M6 thread



2 PCS

6.1805.010

FEP tubing / M6 / 13 cm



2 PCS

6.1805.080

FEP tubing / M6 / 25 cm



2 PCS

6.1805.100

FEP tubing / M6 / 40 cm



1 PCS

6.1829.010

FEP aspiration tube M6, 0.25 m



4 PCS

6.2726.090 拡散防止バルブ

拡散防止用ビュレットチップ6.1541.XXX、 6.1543.130、
6.1543.140、 6.1543.150、 6.1543.160、 6.1543.200

