



Ionizer for 864 Robotic Balance Sample Processor, 230V

2.136.0306

Ionizer (HAUG) for 864 Robotic Balance Sample Processor with power supply for 230V.

Phạm vi giao hàng 2.136.0306

Qt.	Order no.	Description
1 PCS	1.136.0306	Ionizer for 864 Robotic Balance Sample Processor, 230V
		Ionizer (HAUG) for 864 Robotic Balance Sample Processor with power supply for 230V.



1 PCS

6.9013.000

Ionisation rod for ionizer

The ionisation rod is mounted with the holder (6.2057.160) to the 864 Robotic Balance Sample Processor.



Phụ kiện tùy chọn

Order no.	Description	
6.2057.160	Holder for ionization rod for 864 Robotic Balance Sample Processor	
