



USB thermal printer Neo's

2.141.0100

Compact printer with USB interface for the Ti-Touch, the Titrino plus and Dosimat plus and the 91X Meter instruments.

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

Phạm vi giao hàng 2.141.0100

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS	1.141.0100	USB thermal printer Neo's
-------	------------	---------------------------

Compact printer with USB interface for the Ti-Touch, the Titrino plus and Dosimat plus and the 91X Meter instruments.

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.



1 PCS

6.2151.120

Cable USB MINI (OTG) - USB B, 1.8 m

For connecting USB instruments.



Phụ kiện tùy chọn

Order no.	Description	
6.9001.050	Thermal paper for Custom Neos-N / Q3X	
	Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X. New: 25-year legibility with correct storage.	
	1 roll	

