



700 Dosino

2.700.0010

Drive for dosing unit, with permanently installed cable and 9-pin plug (DB9).

This product is no longer available!

Phạm vi giao hàng 2.700.0010

Qt.	Order no.	Description
1 PCS	1.700.0010	700 Dosino
		Drive for Dosing Unit, with fixed cable and 9-pin plug (DB9).



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.



1 PCS

6.2134.200

Adapter DB9 St. /Mini Din8 Bu



1 PCS

6.2148.040

Adapter for Dosino 700 (series 23 or higher)



Phụ kiện tùy chọn

Order no.	Description	
6.2134.020	Adapter DB9 female/mini-DIN8 male DB9 socket and mini-DIN 8 plug. Connection 700.0010 Dosino to 736, 751, 758, 799 Titrinos, 756, 831 KF Coulometer, 711 Liquino.	