



700 Dosino

2.700.0020

The 700 Dosino is a drive for dosing units. With a permanently installed cable and Mini DIN plug.

This product is no longer available!

Phạm vi giao hàng 2.700.0020

Qt.	Order no.	Description
1 PCS	1.700.0020	700 Dosino
		The 700 Dosino is a Dosing Unit drive. With fixed cable and mini-DIN plug.



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.



1 PCS

6.2148.040

Adapter for Dosino 700 (series 23 or higher)



Phụ kiện tùy chọn

Order no.	Description	
6.2134.010	Adapter DB9 male / Mini-DIN8 female DB9 plug and mini-DIN 8 socket. E.g. for the connection of 700.0020 Dosinos with mini-DIN 8 plugs to 726, 796 Titroprocessors.	
