



703 Ti Stand with stirrer and pump

2.703.0010

Compact titration stand with magnetic stirrer and built-in pump for addition of solvents and for aspirating off titration vessel contents.

Phạm vi giao hàng 2.703.0010

Qt.	Order no.	Description
1 PCS	1.703.0010	703 Ti Stand with stirrer and pump
		Compact titration stand with magnetic stirrer and built-in pump for addition of solvents and for aspirating off titration vessel contents.



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1446.060

Stopper / B-14/15 / M10



1 PCS

6.1543.110

Dosing tip / M8 thread

Transfer tip with M8 thread



2 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1608.023

Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



1 PCS

6.1608.030

Round glass bottle / 1000 mL / GL 45



1 PCS

6.1609.010

Adsorber tube with tubing nipple

Lid with tubing nipple.



2 PCS

6.1801.000

PVC tubing / 4 mm / 7 mm / 1 m



1 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1818.000

Aspiration tubing

For 681 Pump Unit and 703 Ti Stand



1 PCS

6.1819.030

Tube

For M8 threaded holes. For 6.1602.1XX bottle attachments and 6.1608.0XX bottles



1 PCS

6.2001.010

Support stand with support rod

For Dosimats and Titrinos



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2108.100

Connection cable for 622/649/703/728, 0.5 m

For supplying power to 622, 649, 728 Stirrers and 703 Ti Stand from control instruments and corresponding Dosimats



1 PCS

6.2739.000

Wrench

For tightening connectors

