

Pump Unit 772 "Aspiration"

2.772.0020

Peristaltic pump for the aspiration of titrated solution and suspensions. The "Aspiration" version contains extensive accessories for aspiration. The pump can be actuated through a Relay Box 731. It is suitable for the conveyance of solvents.

Phạm vi giao hàng 2.772.0020

Qt.	Order no.	Description	
1 PCS	6.1446.160	Thermometer stopper, B-14/15, oblique insertion To have the inserts of the 6.1458.010 titration head in a sloping position.	
1 PCS	6.1543.170	Aspiration tip / M8 thread Aspiration tip for the Sample Processors.	
1 PCS	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	

1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm

With protection against kinks.



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



1 PCS

6.2621.140

Hexagon key 2.5 mm



Phụ kiện tùy chọn

Order no.	Description	
6.1826.160	Viton pump tubing	
	For peristaltic pumps.	