

Pump Unit 772 "Rinse"

2.772.0030

Peristaltic pump for rinsing sample beakers and electrodes. The "Rinse" version contains extensive accessories for rinsing. The pump can be actuated through a Relay Box 731. It is suitable for the conveyance of solvents.

Phạm vi giao hàng 2.772.0030

Qt.	Order no.	Description	
1 PCS	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	
3 PCS	6.1805.060	FEP tubing / M6 / 60 cm With light and kink protection	
1 PCS	6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	

2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



Phụ kiện tùy chọn

Order no.	Description	
6.1826.160	Viton pump tubing	
	For peristaltic pumps.	