



772 Pump Unit

2.772.0110

Peristaltic pump for direct connection to Sample Processors.

Phạm vi giao hàng 2.772.0110

Qt.	Order no.	Description
1 PCS	1.772.0110	772 Pump Unit Peristaltic pump for direct connection to Sample Processors.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



Phụ kiện tùy chọn

Order no.	Description	
6.1826.160	Viton pump tubing For peristaltic pumps.	
6.2160.010	Adapter cable / 2 m / M8 socket / 2 x plug B Adapter cable / 2 m / M8 socket / 2 x plug B. For 772 Pump Unit.	