



786 Swing Head with titration swing arm, left- or right-swinging

2.786.0030

Additional equipment for Robotic Sample Processors XL for direct titration in sample vessels of 75 mL and more.

Phạm vi giao hàng 2.786.0030

Qt.	Order no.	Description
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.050

Robotic arm with titration head for 786 Swing Head,
left or right swinging

Titration head for Swing Head



1 PCS

6.1808.040

Thread adapter M6 outer / M8 inner

Outer thread M6, inner thread M8.



1 PCS

6.1909.030

Stirring propeller / 95 mm

Propeller stirrer with PP shaft, PVDF tip. Length from the SGJ: 104 mm. For
KF applications with 722, 802 Rod Stirrer



2 PCS

6.2042.030

Buret tip clip (rod stirrer)

Clip for proper alignment of buret tips, dosing or aspiration tubings using
a rod stirrer in the titration beaker of a sample changer



