



# 800 Dosino

2.800.0020

Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).

## Phạm vi giao hàng 2.800.0020

| Qt.   | Order no.  | Description                      |
|-------|------------|----------------------------------|
| 1 PCS | 1.800.0040 | 800 Dosino (cable length 0.65 m) |

Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.

