



804 Ti Stand without stand rod

2.804.0010

Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.

Phạm vi giao hàng 2.804.0010

Qt.	Order no.	Description
1 PCS	1.804.0010	804 Ti Stand without stand rod
		Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.



Phụ kiện tùy chọn

Order no.	Description	
6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrando	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrando, Mini-DIN socket - mini-DIN plug	