


824 Easy Sample Changer

2.824.0020

The economic sample changer for routine use. Easy Sample Changer Micro: sample changer for titrations with smaller sample volumes. Up to 24 samples with a maximum volume of 75 mL can be processed in a single run. With the KFT equipment the Easy Sample Changer can be converted to a Karl Fischer sample changer.

Phạm vi giao hàng 2.824.0020

Qt.	Order no.	Description	
1 PCS	6.2034.020	Mounting bracket for Magnetic Stirrer 741 To attach the magnetic stirrer to the «824 Easy Sample Changer	
1 PCS	6.2752.010	Plug cover Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor	
1 PCS	1.741.0010	741 Magnetic Stirrer Magnetic stirrer for sample changers. The magnetic stirrer is located beneath the sample rack.	

24 PCS

6.1432.210

Sample beaker, clear glass, 75 mL, 1 piece

Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



1 PCS

6.1453.250

Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece



1 PCS

6.1458.020

Titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



1 PCS

6.1805.110

FEP tubing / M6 / 80 cm

With light and kink protection



24 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2041.340

Sample rack, 24 x 75 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 24 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2621.140

Hexagon key 2.5 mm




1 PCS

6.2751.080

Safety shield

Safety shield for single-tower Sample Processors without 786 Swing Head.



Phụ kiện tùy chọn

Order no.	Description
6.1453.220	Sample beaker, plastic (PP), 200 mL, 1 piece Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece

