

pH Lab Irda 120 V US

2.827.0015

Phạm vi giao hàng 2.827.0015

Qt.	Order no.	Description	
1 PCS	6.2161.020	Plug-in power supply unit US 120V/6V DC Power supply for 827 pH lab	
1 PCS	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.	
1 PCS	6.2016.050	Support rod / 300 mm	

1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

