

# pH Lab Irda 240 V AUS

2.827.0017

## Phạm vi giao hàng 2.827.0017

| Qt.   | Order no.  | Description                                                             |                                                                                       |
|-------|------------|-------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 1 PCS | 6.2161.040 | Plug-in power supply unit AUS 240V/6V DC<br>Power supply for 827 pH lab |    |
| 1 PCS | 6.2013.010 | Clamping ring<br>For support rods with a diameter of 10 mm.             |   |
| 1 PCS | 6.2016.050 | Support rod / 300 mm                                                    |  |

1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



---

1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

