

pH Lab Irda 230 V UK

2.827.0019

Phạm vi giao hàng 2.827.0019

Qt.	Order no.	Description
1 PCS	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.
		
1 PCS	6.2016.050	Support rod / 300 mm
		
1 PCS	6.2021.020	Electrode holder Electrode holder for 4 electrodes and 2 buret tips
		

1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2161.010

Plug-in power supply unit EU 230V /6V DC

Power supply for 827 pH lab

