



827 pH lab; 240 V, AUS with Primatrode

2.827.0117

Laboratory pH meter with IrDa interface and 6.0228.010 Primatrode

Phạm vi giao hàng 2.827.0117

Qt.	Order no.	Description
1 PCS	6.2161.040	Plug-in power supply unit AUS 240V/6V DC
		Power supply for 827 pH lab



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



Phụ kiện tùy chọn

Order no.	Description	
6.2307.230	Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	pHit kit Maintenance kit for pH electrodes The kit contains: • 50 mL cleaning solution • 50 mL 3M KCl solution • 50 mL storage solution • 2 Storage vessels • Instructions for use	