



# 827 pH lab; 230 V, EU with Unitrode

2.827.0214

Laboratory pH meter with IrDa interface and 6.0258.600 Unitrode

Phạm vi giao hàng 2.827.0214

| Qt.   | Order no.  | Description          |
|-------|------------|----------------------|
| 1 PCS | 6.0258.600 | Unitrode with Pt1000 |

Combined pH electrode with integrated Pt1000 temperature sensor for pH measurements/titrations. This electrode is particularly suitable:

- for pH measurements and titrations in difficult, viscous, or alkaline samples
- at elevated temperatures
- for long-term measurements



The fixed ground-joint diaphragm is insensitive to contamination.

Reference electrolyte:  $c(\text{KCl}) = 3 \text{ mol/L}$ , storage in storage solution.

Alternatively: reference electrolyte for measurements at  $T > 80^\circ\text{C}$ : Idrolyte, storage in Idrolyte.

1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



---

1 PCS

6.2016.050

Support rod / 300 mm



---

1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



---

1 PCS

6.2104.600

Electrode cable for plug in head U/plug F, 2x2 mm B,  
1m

For connecting electrodes with Metrohm plug-in head U to Metrohm  
instruments (socket F).



1 PCS

6.2161.010

Plug-in power supply unit EU 230V /6V DC

Power supply for 827 pH lab



---

1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



## Phụ kiện tùy chọn

| Order no.  | Description                                                                                                                                                                                                                                                                      |                                                                                     |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2307.230 | Buffer solutions pH 4, 7 and 9<br><br>Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL                                                                                                                                  |  |
| 6.2325.000 | pHit kit<br><br>Maintenance kit for pH electrodes<br><br>The kit contains: <ul style="list-style-type: none"><li>• 50 mL cleaning solution</li><li>• 50 mL 3M KCl solution</li><li>• 50 mL storage solution</li><li>• 2 Storage vessels</li><li>• Instructions for use</li></ul> |  |