



833 Advanced IC Liquid Handling Suppressor Unit

2.833.0020

Chemical «packed bed anion suppressor. Patented MSM Metrohm Suppressor Module as stand-alone instrument or for retrofitting HPLC or IC systems. Solvent- and pressure-resistant. With two-channel pump for regeneration and rinsing solution.

Phạm vi giao hàng 2.833.0020

Qt.	Order no.	Description
1 PCS	6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m
		For all IC instruments



2 PCS

6.2128.180

Cable 762/830 - 752/753/754/833

Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit



2 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.




2 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



2 PCS



6.2821.120

Inline filter 2 µm

Inline filter 2 µm



Phụ kiện tùy chọn

Order no.	Description	
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
6.1826.060	Pump tubing (orange/yellow) For Compact IC , Personal IC , Basic IC	
6.2755.000	Tubing cartridge For the peristaltic pumps of IC instruments. Contact pressure adjustable.	
6.2821.130	Spare filter for inline filter Spare filter for inline filter.	

6.2832.000

MSM Rotor A

Suppressor rotor for all IC instruments with MSM (Metrohm Suppressor Module)



6.2832.010

Connecting piece for MSM and MSM-LC

With inlets and outlets. For all IC instruments with MSM (Metrohm Suppressor Module).

For replacement in Professional IC instruments of generations 850, 881, 882, and 883 (manufactured before 2017).

