



833 IC Liquid Handling Ultra Filtration Unit

2.833.0050

For inline ultrafiltration. Guarantees the reliable removal of particles from the sample. 2-channel peristaltic pump and separate ultrafiltration cell, for example for installation in the 820 IC Separation Center.

Phạm vi giao hàng 2.833.0050

Qt.	Order no.	Description
1 PCS	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC.



1 PCS

6.1803.050

PTFE capillary 0.5 mm i.d. / 20 cm

Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor



2 PCS

6.1803.070

PTFE capillary 0.97 mm i.d. / 40 cm

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



2 PCS

6.1803.080

PTFE capillary 0.97 mm i.d. / 1 m

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



1 PCS

6.1803.100

PTFE capillary, 0.5 mm ID / 40 cm

Capillary for filtration and dialysis



1 PCS

6.2128.180

Cable 762/830 - 752/753/754/833

Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces

Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



1 PCS

6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.030

Coupling olive/UNF 10/32, 4x

Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 PCS


6.2831.010

Plastic tweezers

For dialysis and filtration membranes



Phụ kiện tùy chọn

Order no.	Description	
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
6.1826.060	Pump tubing (orange/yellow) For Compact IC , Personal IC , Basic IC	
6.2755.000	Tubing cartridge For the peristaltic pumps of IC instruments. Contact pressure adjustable.	