



843 Pump Station (membrane) - rinse /aspirate for Sample Processors

2.843.0050

The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly from the Sample Processor via pump commands. The Rinse / Aspirate version is provided with all the accessories needed for automatically emptying the titration beaker and rinsing the titration equipment. The accessories are intended for use with the 814 and 815 Sample Processors.

Phạm vi giao hàng 2.843.0050

Qt.	Order no.	Description
1 PCS	1.843.0020	843 Pump Station (membrane)
The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.		



1 PCS

6.1543.170

Aspiration tip / M8 thread

Aspiration tip for the Sample Processors.



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1621.100

Cap for 10-L container

Cap for 10-L container. For easy and safe transportation of filled containers.



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS 6.1805.060 FEP tubing / M6 / 60 cm
With light and kink protection



3 PCS 6.1805.420 FEP tubing / M6 / 48 cm
With light and kink protection



1 PCS 6.1805.510 PTFE tubing / M8 / 60 cm
With protection against kinks.



1 PCS 6.1808.170 Distributor for Sample Processors for rinsing and
aspirating

Distributor for installation on Sample Processors for rinsing and aspirating
when using externally connected pumps.



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



1 PCS

6.2141.300

Remote cable 843 Pump Station to Sample Processor

Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



Phụ kiện tùy chọn

Order no.	Description	
6.2740.030	Rinsing nozzle The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.	
