



891 Professional Analog Out

2.891.0010

The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.

Phạm vi giao hàng 2.891.0010

Qt.	Order no.	Description
1 PCS	1.891.0010	891 Professional Analog Out

The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.



1 PCS

6.2128.200

Cable 891-Analog

Analog cable to connect to external AD converters.



1 PCS

6.2156.070

Cable Extension Module - Professional IC, 1 m

Cable connecting the extension module with the Professional IC instrument

