



Generator electrode without diaphragm

6.00349.100

Generator electrode without diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well-suited for most applications where the sample is highly soluble in alcohol.

This electrode has a cable connection angled upwards.

Phụ kiện tùy chọn

Order no.	Description
6.2104.120	Electrode cable / 1 m / H With plug H. For connecting generator electrode - KF Coulometer.



6.2104.620

Electrode cable / 0.65 m / H

With plug H, for the generator electrode - coulometer connection

