



Measuring head (scTRACE Gold) for the 946 Portable VA Analyzer measuring stand (for older instrument versions)

6.01256.000

Measuring head with built-in stirrer for measurements with the scTRACE Gold.

Phụ kiện tùy chọn

Order no.	Description
6.01204.000	Stirrer tip
	Stirrer tip for 946 Portable VA Analyzer.



6.1258.000

scTRACE Gold (set with 4 pieces)


Gold micro-wire electrode on thick film basis for the determination of arsenic in waters. The scTRACE Gold is a combined voltammetric sensor for voltammetry. Gold micro-wire working electrode, Ag/AgCl reference electrode and carbon auxiliary electrode are located on a sensor. Simple in its handling and practically maintenance-free, the scTRACE Gold can for example be used to determine arsenic in water in the trace range. Together with the 6.1241.080 electrode shaft, it can be used with any Metrohm voltammetry measuring stand.



6.01408.000

O-ring 46 x 2.5 mm

For the 946 Portable VA Analyzer measuring stand. Set of 2 pieces.



6.02621.000

4 mm wrench

