



## Measuring head (SPE) for the 946 Portable VA Analyzer measuring stand

6.01256.020

Measuring head with integrated stirrer for measurements with *screen-printed electrodes* (SPE).

### Phụ kiện tùy chọn


| Order no.   | Description |
|-------------|-------------|
| 6.01204.000 | Stirrer tip |

Stirrer tip for 946 Portable VA Analyzer.



6.01408.000 O-ring 46 x 2.5 mm

For the 946 Portable VA Analyzer measuring stand. Set of 2 pieces.



---

6.1241.200 Nut for SPE electrode shaft

Nut for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020.



---

6.1241.210 Supporting ring for SPE electrode shaft

Supporting ring made of PTFE for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020.



---

6.1244.060 Silicone seal for SPE electrode shaft (10 pieces)

Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.



6.02621.000

4 mm wrench

