



Titration head 3xNS14 / 4x6.4 mm (P&P)

6.01403.030

Titration head with 3 NS14 and 4x 6.4 mm openings for using 120 mL beakers in Pick&Place modules of the OMNIS Sample Robot Pick&Place with two retracted dosing tips. Prepared for being equipped with rinsing and aspiration tubing.

Phạm vi giao hàng 6.01403.030

Qt.	Order no.	Description
2 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



2 PCS

6.1446.000

SGJ stopper / B-14/(15)



2 PCS

6.1808.000

Connection bushing M6



2 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Phụ kiện tùy chọn

Order no.	Description	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	