



OMNIS titration vessel lid for 601406XXX titration vessels for Polytron

6.01405.040

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator with Polytron.

Phạm vi giao hàng 6.01405.040

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments



Phụ kiện tùy chọn

Order no.	Description	
6.02709.030	Screw nipple Plastic screw nipple, suitable for 6014050X0 titration vessel lids.	
6.02709.020	Septum stopper Plastic septum stopper, suitable for 6014050X0 titration vessel lids.	
6.02709.040	Stopper Stopper for titration vessel lids 6014050X0.	