



Titration head adapter

6.01462.010

This adapter enables the utilization of titration heads of the Sample Processor generation (814 /815/855) on sample changers with screw-on titration heads.

Phụ kiện tùy chọn

Order no.	Description
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings. With 6xSGJ 14 and 3x SGJ 9 openings.



6.1458.020

Titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



6.1458.030

KF titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



6.1458.040

Titration Head, 3 x SGJ 14

Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.



6.1458.070

Titration head / 2 x SGJ 14

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.

