



Cylinder unit OMNIS special 20 mL, without accessories

6.01508.220

Intelligent cylinder unit 20 mL for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions

Without dosing tubing and antidiffusion tip.

Phạm vi giao hàng 6.01508.220

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

