



OMNIS sample rack, 9 x 150 mL, (PP)

6.02041.050

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 9 sample beakers. The following sample beakers can be used: 150 mL beaker (tall shape and without spout).

Plastic: Polypropylene (PP)

Phạm vi giao hàng 6.02041.050

Qt.	Order no.	Description
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Phụ kiện tùy chọn

Order no.	Description
6.01404.050	Beaker adapter for 150 mL sample beaker Beaker adapter for centering the 150 mL sample beaker in the slide of the Pick&Place modules. 

