




# OMNIS sample rack 77 x 8 mm

6.02041.080

OMNIS sample rack for OMNIS Sample Robot – NIR suitable for 77 disposable vials with an optical path length of 8 mm for measurements in transmission (6.7402.240).

Phạm vi giao hàng 6.02041.080

| Qt.   | Order no.   | Description                |
|-------|-------------|----------------------------|
| 1 PCS | 6.02042.000 | Holder for rack label (3x) |

The holder for rack labels permits durable labeling of the rack with an unambiguous name.

