




KF Dis-Cover lid for OMNIS 75 mL sample beaker, 25 pieces

6.02710.000

Lid with seals for volumetric Karl Fischer Titration in 6.01402.000 and 6.01402.003 sample beakers in the OMNIS Sample Robot Pick&Place system.

Phạm vi giao hàng 6.02710.000

Qt.	Order no.	Description
1 PCS	6.05700.160	Consumable Kit seals for the 6.02710.000 KF Dis-Cover lid

Replacement seals for 6.02710.000 lid, 25 pieces.

