



Ultrafiltration cell 2

6.02729.110

For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.

Phụ kiện tùy chọn

Order no.	Description
6.2714.020	Filtration membrane reg. cellulose, 0.2 µm, 50 pieces
	Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



6.2714.030 Dialysis membrane (polyamide)
For inline dialysis and filtration



6.2744.000 Pressure screw PVDF
With UNF 10/32 connection. For IC instruments, inline dialysis
(connection of PTFE capillaries)



6.2831.010 Plastic tweezers
For dialysis and filtration membranes



6.02057.120 Ultrafiltration cell 2 holder, stand-alone
Holder for the 6.02729.110 ultrafiltration cell 2



6.02057.030 Ultrafiltration cell 2 holder

Holder for mounting the 6.02729.110 ultrafiltration cell 2 on the 858 Professional Sample Processor / 919 IC Autosampler plus.



6.2714.050 Filtration membrane reg. cellulose, 0.45 µm, 50 pieces

Nominal pore size: 0.45 µm. Set of 50 pieces. For inline filtration.

