



Adapter Luer inner / M6 outer

6.02744.050

Adapter for syringe for pre-rinsing the aspiration filter (6.2821.090).

The adapter can be screwed onto the aspiration filter. A syringe can be used to rinse the aspiration filter with ultrapure water or eluent.

Phụ kiện tùy chọn

Order no.	Description
6.2821.090	Aspiration filter
	Pore size: 27 µm. Set of 5 pieces. For aspiration tubing (6.1834.000) and filter tubes (6.1821.040) and (6.1821.050).



6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA

