



Cylinder unit OMNIS 10 mL

6.03001.210

Intelligent 10 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.

Phạm vi giao hàng 6.03001.210

Qt.	Order no.	Description
1 PCS	6.01503.210	OMNIS 10 mL cylinder unit without accessories

Intelligent 10 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



1 PCS

6.1446.030

Ball stopper / B-14/(15)

For the movable holding of buret tips in titration vessels



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



2 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.2008.030

Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



1 PCS


6.2739.000

Wrench

For tightening connectors



Phụ kiện tùy chọn

Order no.	Description	
6.05700.140	Consumable Kit, cylinder unit OMNIS	
	FEP tubing connections M6 (40 cm) and antidiﬀusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.	
6.05700.421	Consumable Kit, cylinder unit OMNIS 10 mL	
	Cylinder unit OMNIS 10 mL complete with cylinder, piston, cylinder base, valve disk, centering tube, and cylinder housing, including programmed data chip.	