



# Cylinder unit OMNIS 50 mL

6.03001.250

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.

Phạm vi giao hàng 6.03001.250

| Qt.   | Order no.   | Description                                   |
|-------|-------------|-----------------------------------------------|
| 1 PCS | 6.01503.250 | OMNIS 50 mL cylinder unit without accessories |

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



1 PCS

6.1446.030

Ball stopper / B-14/(15)

For the movable holding of buret tips in titration vessels



---

2 PCS

6.1446.040

Threaded stopper / M6



---

1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



---

2 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.2008.030

Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



---

1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



---

1 PCS

6.2739.000

Wrench

For tightening connectors



# Phụ kiện tùy chọn

| Order no.   | Description                                                                                                                                                |                                                                                     |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.05700.140 | Consumable Kit, cylinder unit OMNIS                                                                                                                        |  |
|             | FEP tubing connections M6 (40 cm) and antidiﬀusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.                                               |                                                                                     |
| 6.05700.441 | Consumable Kit, cylinder unit OMNIS 50 mL                                                                                                                  |  |
|             | Cylinder unit OMNIS 50 mL complete with cylinder, piston, cylinder base, valve disk, centering tube, and cylinder housing, including programmed data chip. |                                                                                     |