



# Double Pt-wire electrode for coulometry

6.0341.100

Indicator electrode used for coulometric Karl Fischer titration.

Phạm vi giao hàng 6.0341.100

| Qt.   | Order no.  | Description    |
|-------|------------|----------------|
| 1 PCS | 6.2008.040 | Storage vessel |

Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

| Order no.  | Description                                                                                                               |                                                                                     |
|------------|---------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2104.020 | Electrode cable / 1 m / F<br><br>For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F). |  |