



Micro Pt Titrode (installation length 17.8 cm)

6.0434.110

Combined platinum ring electrode with a pH glass membrane as reference electrode, installation length = 17.8 cm, narrowed to a diameter of 6.4 mm.

This maintenance-free electrode is well suited for redox titrations when the pH value remains constant of, e.g.:

- iodometry
- chromatometry
- cerimetry
- permanganometry

This electrode is stored in deion. water.

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.2104.020

Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).

