



## Combined iPt ring electrode

6.0451.300

Intelligent, combined platinum ring electrode with integrated memory chip for storing sensor data.

This electrode is suitable for redox titrations where the pH value varies, e.g.:

- Oxygen content according to Winkler
- Determination of hydrogen peroxide with  $\text{KMnO}_4$
- Diazotization titrations

$c(\text{KCl}) = 3 \text{ mol/L}$  is used as reference electrolyte and for storage.

iTodes can be used on Titrand, Ti-Touch or 913/914 meters.

Phạm vi giao hàng 6.0451.300

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



---

1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



## Phụ kiện tùy chọn

| Order no.  | Description                                                                                                                                |                                                                                     |
|------------|--------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 2.854.0010 | 854 iConnect<br><br>The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm). |  |
| 6.2308.020 | Electrolyte 3 mol/L KCl (250 mL)<br><br>Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)           |  |